

Complete Cair™ Boot

Indications for Use:

Promotes neutral alignment of the foot and relieves pressure on the heel.

Product Features:



- **Non-slip sole:** Allows short standing and walking.
- **Foam padded insole:** For optimal comfort.
- **Outer shell:** Antimicrobial, wipedown fabric for easy cleaning.
- **Internal low-profile air bladder:** Allows for dynamic contouring and support around back of the leg and shin.
- **Patented Contouring pad:** Redistributes skin pressure on the leg and ankle.
- **Closure Straps:** These flexible straps stretch to comfortably secure the leg.

Instructions for use:



- 1 Remove boot from bag & open straps.

Place the foot into the boot with the heel centered over the opening. The contouring pad will gently cradle the Achilles tendon proximal to the ankle.



- 2 Secure the flap over the shin, starting with the strap closest to the knee. Ensure the flaps are secure, but not too tight. Next, attach the blue talus strap, followed by the strap across the top of the foot.



- 3 Pass the free ends of the black dorsiflexion-assist strap through the D-rings at the top of the boot and attach them onto themselves. Adjust to achieve the desired tension. Tighten one side as to address inversion/eversion. Straps can be omitted or removed.

Complete Cair™ Boot

Product Sizing:

Universal Fit

Single Patient Use:

Product can be wiped down using hospital approved disinfectant.

Additional Considerations:

The Complete Cair Boot should be removed at least once per shift or more as dictated by hospital protocol to inspect the foot and lower leg.

Product Warranty:

 Limited 90-day warranty against manufacture defects.

Product Details:

- This product is not made with natural rubber latex (NRL) or di-(2-ethylhexyl)phthalate (DEHP).



Scan here for full
in-service video

