

Indications for Use:

Promotes neutral alignment of the foot and relieves pressure on the heel.

Features + Benefits:

- **Non-slip sole:** Allows for stand/pivot. Not made for walking.
- **Foam padded insole:** For optimal comfort.
- **Detachable outer shell:** Helps prevent foot drop.
- **Antimicrobial cover:** Wipe down fabric for easy cleaning.
- **Internal low-profile air bladder:** Allows for dynamic contouring and support around back of the leg and shin.
- **Closure straps:** Flexible straps comfortably secure the leg.
- **Gate feature:** Easy access to foot without removing the boot.
- **Detachable anti-rotation wedge:** Prevents internal/external rotation of the foot.



Instructions for Use:

- 1** Remove boot from bag. Place the foot into the boot with the heel centered over the opening. The contouring pad will gently cradle the Achilles tendon proximal to the ankle.
- 2** Close the adjustable straps, beginning with the one closest to the knee and working down towards the foot. The straps are elastic to accommodate the shape of the leg. Make sure they are secure but not too tight. You should be able to slide a finger between the boot and leg to confirm appropriate fit.
- 3** Attach both sides of foot-drop strap with foot in neutral alignment. Make sure the anti-rotation wedge is secured to desired angle.

