


Appendix A - Accessory for NeuroFlex® Restorative Knee

Knee Separator (Product code: 20725SP)

Helps to gradually separate the knee joint before transitioning to the NeuroFlex® Restorative Knee



1. Follow your clinicians' recommendation on the number of foam blocks (1, 2, or none) that should be used.
2. If using foam block(s), open the flap at the back (blue) side of the Knee Separator and insert the block(s). Otherwise, save them for future use.
3. Open the three blue straps.
4. Fold the Knee Separator in half with the black side on the outside and all the straps on the inside of the leg.
5. Apply the folded Knee Separator to the inside of the knee with the black material against the patient's leg, centering the fold against the knee crease.
6. Wrap the straps around the leg, attaching them to the blue fabric so the Knee Separator is firmly but comfortably secured to the patient's leg, as shown in the lower picture.
7. Adjust the Knee Separator to gently increase the knee's range of motion and to abduct the legs/hips if recommended by your clinician.
8. Add extra foam padding into the pocket when recommended by your clinician to continue the improvement.
9. Continue to use the Knee Separator, if recommended by your clinician, after transitioning to the NeuroFlex® Restorative Knee during the period that the latter is not being worn to prevent skin-on-skin contact and to maintain range of motion.