

# Fluidized AFO

## Prevent - Protect - Restore



**Fluidized Technology**

**Thermal Adaptive Technology**

- Safely and comfortably addresses a non-ambulatory solution for contractures of the foot and ankle, foot drop, plantar flexion and conditions resulting in various types of neuropathy, inflammation or nerve injury.
- Patented fluidized cushioning creates the perfect blend of support and contouring capability to help achieve the best possible therapeutic positioning.
- Helps to maintain neutral alignment of the foot in a non-ambulatory patient.
- Fluidized sections to help redistribute force include lower calf and Achilles area, and the counterforce Talus Strap.
- Suspends the heel to prevent pressure and skin breakdown.
- Anti-skid pad for safe standing and pivot transfers (not designed for ambulation).
- High temperature plastic is molded to help float the heel, and is easily remolded if needed to accommodate moderate to severe lost range of motion.

AFO Boot Size	Foot Size
Small	Up to 8"
Medium	Up to 9.5"
Large	Up to 11"
X-Large	Up to 13"

- 3 dimensional Fluidized Medium micro contours at the Achilles area and reduces friction and shear
- Soft, durable fabric next to patient's skin has Thermal Adaptive Technology to help keep foot cool and dry
- Rotator Bar helps maintain alignment of lower leg, helps to correct hip rotation, and is completely covered with padded material
- Adjustable hook and loop enclosures allow caregiver to customize for a more perfect fit
- Adjustable straps maximize dorsiflexion to help prevent foot drop
- Available in small, medium, large, and x-large - custom sizes available upon request



Prevent • Protect • Restore

332 East Broadway • Brandenburg, KY 40108  
(800) 793-5544 • 270-422-5454 • Fax: (270) 422-5453

[www.restorativemedical.com](http://www.restorativemedical.com)



RM- AFO Sales Sheet v1.6 - 10/31/2016