

ELBOW / SHOULDER



Measuring Guide

NeuroFlex® Restorative Elbow Item Code	NeuroFlex® Restorative HyperExtension Elbow Item Code	Size	BTE	ETF	Biceps Cir.	Forearm Cir.
20724C	20724H	Child	3 ½"	4"	5 ½-7 ½"	3 ½-6 ½"
20724CF	20724CFH	Child/Full	3 ½"	4"	6 ½-8 ½"	5 ½- 7 ½"
20724YS	20724YSH	Youth/Slim	4"	5"	5 ½-7 ½"	3 ½-6 ½"
20724Y	20724YH	Youth	4"	5"	6 ½-8 ½"	5 ½- 7 ½"
20724XS	20724XSH	Ex-Small	5"	5"	7 ½-9 ½"	6 ½-8 ½"
20724XSF	20724XSFH	Ex-Small/Full	5"	5"	8-12"	7-10"
20724SS	20724SSH	Small/Slim	5"	6"	7 ½-9 ½"	6 ½-8 ½"
20724S	20724SH	Small	5"	6"	8-12"	7-10"
20724SF	20724SFH	Small/Full	5"	6"	8-13"	7 ½-11"
20724MS	20724MSH	Medium/Slim	6"	6"	8-12"	7-10"
20724M	20724MH	Medium	6"	6"	8-13"	7 ½-11"
20724MF	20724MFH	Medium/Full	6"	6"	10-15"	8-12"
20724LS	20724LSH	Large/Slim	6 ½"	6 ½"	8-13"	7 ½-11"
20724L	20724LH	Large	6 ½"	6 ½"	10-15"	8-12"
20724LF	20724LFH	Large/Full	6 ½"	6 ½"	13-17"	10-15"
20724Z	20724ZH	Custom				