



Measuring Guide

RM Pediatric HyperHand Item Code	NeuroFlex® Restorative ThumbEase (Pediatric) Item Code	Size	Width of MCP	Wrist to Finger	Wrist to Forearm
20717IR	N/A	Infant-Right	1 ¾" to 2"	Up To 2 ½"	2 ¼"
20717IL	N/A	Infant-Left	1 ¾" to 2"	Up To 2 ½"	2 ¼"
20717TJR	N/A	Toddler Jr.-Right	1 ¾" to 2"	Up To 3"	2 ¾"
20717TJL	N/A	Toddler Jr.-Left	1 ¾" to 2"	Up To 3"	2 ¾"
20717TR	20740TR	Toddler-Right	2" to 2 ¼"	Up To 3 ½"	3 ½"
20717TL	20740TL	Toddler-Left	2" to 2 ¼"	Up To 3 ½"	3 ½"
20717CR	20740CR	Child-Right	2 ¼" to 2 ½"	Up To 4"	3 ¾"
20717CL	20740CL	Child-Left	2 ¼" to 2 ½"	Up To 4"	3 ¾"
20717YR	20740YR	Youth-Right	2 ¾" to 3"	Up To 5"	4"
20717YL	20740YL	Youth-Left	2 ¾" to 3"	Up To 5"	4"