

# RESTORATIVE MEDICAL PRODUCT SUMMARY

All Restorative Medical splints can be adjusted as the patient improves by bending, heat molding, or adjusting hinges.

## HAND



**RM Restorative Hand**  
**Product # 20736**  
**L-Code L3807/L3809**  
Comfortably addresses wrist/hand/finger lost range and tension. Helps to relax muscle tone and spasticity.



**RM BendEase Hand - Universal**  
**Product # 20715**  
**L-Code L3807/L3809**  
Moldable aluminum base, used for moderate lost range (universal size fits most adults).



**RM Cock-Up Wrist Splint**  
**Product # 20739**  
**L-Code L3807/L3809**  
Provides moderate, comfortable support of hand and wrist post-injury or surgery, with 30 degree wrist cock-up.



**RM Cock-Up Hand Splint with Finger Extension Roll**  
**Product # 20739roll**  
**L-Code L3807/L3809**  
30 degree wrist cock-up with removable large roll when finger extension is desired.



**NeuroFlex® Restorative Flex Hand**  
**Product # 20721**  
**L-Code L3807/L3809**  
Addresses tightest closed fist to moderate finger extension with graduated size cones, reverse cones or cylinders.



**NeuroFlex® Restorative ThumbEase**  
**Product # 20737**  
**L-Code L3807/L3809**  
For severely adducted thumbs and thenar eminence collapse.



**RM BendEase Hand - Sized**  
**Product # 20715**  
**L-Code L3807/L3809**  
Bend-by-hand moldable aluminum base, used for moderate lost range.



**RM Prosperity Hand**  
**Product # 20732**  
**L-Code L3807/L3809**  
Slit between the 3rd & 4th digits allows for varying degrees of range of motion in fingers.



**RM Restorative Wrist**  
**Product # RM7001**  
**L-Code: L3915/3916**  
Provides stability and support for the wrist. Adjustable wrist flexion and extension. Removable finger separators.

**All Restorative Medical products are created specifically to:**

### PREVENT

Long-term complications, further injury to joints and surrounding tissue, as well as complications resulting from poor alignment and restricted movement

### PROTECT

Skin from excessive pressure and shearing

### RESTORE

Lost range of motion and function

## FOOT/ANKLE



**NeuroFlex® Restorative FootEase**  
**Product # 20742**  
**L-Code L4396/4397**  
Addresses mild to moderate plantar flexion contractures, extensor tone and plantar fasciitis.



**RM Safe Boot II**  
**Product # 20744II**  
**L-Code L4396/L4397**  
Addresses plantar flexion, flaccid foot drop, plantar fasciitis and mild abnormal muscle tone.



**RM Safe Boot I**  
**Product # 20744I**  
**L-Code L4398**  
Addresses flaccid foot drop and suspends the heel.



**RM Fluidized Boot**  
**Product # 20743**  
**L4396/4397**  
Helps to maintain neutral alignment of the foot and to relieve pressure on the heel using fluidized technology.



**Flo-Form® Complete Cair Boot**  
**Product # RM20747U**  
**L-Code L4396/L4397**  
Helps maintain neutral alignment of the foot and to relieve pressure on the heel.



**NeuroFlex® Restorative Flex Boot**  
**Product # 20734**  
**L-Code L4396/4397**  
Addresses mild to severe lost range and abnormal muscle tone.



**EZ Stride Anterior Carbon Fiber AFO**  
**Product # CFA**  
**L-Code L1932**  
Enhances normal gait for persons with flaccid foot drop. The Lateral strut design spirals to anterior aspect, providing support for patients who pronate (valgus or eversion). Tibial pad available.



**EZ Stride Posterior Carbon Fiber AFO**  
**Product # CFP**  
**L-Code L1951**  
Enhances normal gait for persons with flaccid foot drop. The lateral strut design spirals to medial aspect. Best for patients who cannot tolerate tibial pressure.

## ELBOW/SHOULDER



**NeuroFlex® Restorative Elbow**  
**Product # 20724**  
**L-Code L3760**  
Addresses mild to severe flexion contractures. Available in styles to treat extensor tone/hyperextended elbows, and for below elbow amputations (BEA).



**Elbow Separator (Accessory for NeuroFlex® Restorative Elbow)**  
**Product # 20724SP**  
Helps to gradually separate the elbow joint before transitioning to the NeuroFlex® Restorative Elbow.



**RM RestAir Elbow**  
**Product # 20724air**  
**Suggested L-Code L3760**  
Addresses elbow flexion contractures with air technology.



**NeuroFlex® Restorative Hyper-Extension Elbow**  
**Product # 20724H**  
**L-Code L3760**  
Helps to improve and/or restore extension lost range of the elbow



**Restorative Shoulder**  
**Product # RM6001**  
**L-Code L3960**  
Provides stability and support for chronic shoulder pain.

# RESTORATIVE MEDICAL PRODUCT SUMMARY

Our comprehensive selection of static splints incorporate dynamic properties and technologies to provide therapeutic positioning that will aid in the management and rehabilitation of contractures. Sizes range from infant to adult.



Prevent • Protect • Restore

restorativemedical.com

BACK/NECK



**RM VertebrEase TLSO**  
**Product # RM20730**  
**L-Code L0456/L0457**  
Helps improve posture for spinal deformities such as: kyphosis, lordosis, scoliosis and more. Padded shoulder straps encourage upright posture.



**RM VertebrEase LSO**  
**Product # 20729**  
**L-Code L0631/L0648**  
Helps immobilize and support the lumbar spine for soft tissue injuries and relieve lower back pain. Shown with optional shoulder straps. Available with pendulous abdomen front.



**NeuroFlex® Restorative Safe Spine**  
**Product # 20728**  
**L-Code L0454/L0455**  
Comfortable spine support for conditions like kyphosis, scoliosis, lordosis. Padded shoulder straps encourage upright posture.



**NeuroFlex® Restorative Cervical Extension, Forward**  
**Product # 20728C**  
**L-Code L0999**  
Fits into our VertebrEase and Safe Spine TLSOs and addresses forward head upper body posture.



PDAC Approved

**Restorative LSO 637**  
**Product # RM20637**  
**L-Code L0637/L0650**  
Helps to improve posture and lumbar spinal stability for herniated disc, stenosis, spondylosis, degenerative disc disease, facet joint syndrome, and post-operative rehabilitation.



PDAC Approved

**Restorative TLSO 457**  
**Product # RM20456**  
**L-Code L0457**  
Helps to improve posture and spinal stability for herniated disc, stenosis, spondylosis, degenerative disc disease, facet joint syndrome, and post-operative rehabilitation.



**Restorative Lumbar Sacral Orthosis (LSO) Stringback**  
**Product # RM10037**  
Provides support, helps manage trunk stabilization, and improves posture and torso alignment



**NeuroFlex® Restorative Cervical HyperExtension Stop**  
**Product # 20728H**  
**L-Code L0999**  
Fits into our TLSO and addresses neck hyper-extension, head leaning.



**NeuroFlex® Restorative Safe Cerv, Left or Right**  
**Product # 20728SCL/R**  
**L-Code L0999**  
Fits into our TLSO and addresses left or right head leaning with or without forward leaning.



**NeuroFlex® Restorative Kentucky Kollar**  
**Product # 20735**  
**L-Code L0113**  
Addresses dystonia, torticollis and wry neck conditions.

KNEE/HIP



**NeuroFlex® Restorative Knee**  
**Product # 20725**  
**L-Code L1831**  
Addresses mild to severe flexion contractures. Available in styles to treat extensor tone/hyperextended knees, and for below knee amputations (BKA).



**Knee Separator (Accessory to NeuroFlex® Restorative Knee)**  
**Product # 20725SP**  
Helps to gradually separate the knee joint before transitioning to the NeuroFlex® Restorative Knee.



**RM RestAir Knee**  
**Product # 20755air**  
**L-Code L1831/L1847**  
Addresses knee flexion contractures with air technology.



**NeuroFlex® Restorative Hyperextension Knee**  
**Product # 20725H**  
**Suggested L-Code L1831**  
Helps to improve and/or restore extension lost range of the knee with effective prolonged low load passive stretch.



**RM RestAir Hip**  
**Product # 20760air**  
**L-Code L1652**  
Comfortably and safely helps to address hip adduction & scissoring.



**Restorative OA Knee**  
**Product # 9021**  
**L-Code: L1843**  
Provides stability and support for the knee.



**Restorative Hinged Knee**  
**Product # 9001**  
**L-Code: L1832/L1833**  
Provides stability and support for chronic knee pains and instability.