

RESTORATIVE MEDICAL PRODUCT SUMMARY

All Restorative Medical splints can be adjusted as the patient improves by bending, heat molding, or adjusting hinges.

HAND



RM Restorative Hand
Product # 20736
L-Code L3807/L3809
Comfortably addresses wrist/hand/finger lost range and tension. Helps to relax muscle tone and spasticity.



RM BendEase Hand - Universal
Product # 20715
L-Code L3807/L3809
Moldable aluminum base, used for moderate lost range (universal size fits most adults).



RM Cock-Up Wrist Splint
Product # 20739
L-Code L3807/L3809
Provides moderate, comfortable support of hand and wrist post-injury or surgery, with 30 degree wrist cock-up.



RM Cock-Up Hand Splint with Finger Extension Roll
Product # 20739roll
L-Code L3807/L3809
30 degree wrist cock-up with removable large roll when finger extension is desired.



NeuroFlex® Restorative Flex Hand
Product # 20721
L-Code L3807/L3809
Addresses tightest closed fist to moderate finger extension with graduated size cones, reverse cones or cylinders.



NeuroFlex® Restorative ThumbEase
Product # 20737
L-Code L3807/L3809
For severely adducted thumbs and thenar eminence collapse.



RM BendEase Hand - Sized
Product # 20715
L-Code L3807/L3809
Bend-by-hand moldable aluminum base, used for moderate lost range.



RM Prosperity Hand
Product # 20732
L-Code L3807/L3809
Slit between the 3rd & 4th digits allows for varying degrees of range of motion in fingers.



RM Restorative Wrist
Product # RM7001
L-Code: L3915/3916
Provides stability and support for the wrist. Adjustable wrist flexion and extension. Removable finger separators.

All Restorative Medical products are created specifically to:

PREVENT

Long-term complications, further injury to joints and surrounding tissue, as well as complications resulting from poor alignment and restricted movement

PROTECT

Skin from excessive pressure and shearing

RESTORE

Lost range of motion and function

FOOT/ANKLE



NeuroFlex® Restorative FootEase
Product # 20742
L-Code L4396/4397
Addresses mild to moderate plantar flexion contractures, extensor tone and plantar fasciitis.



RM Safe Boot II
Product # 20744II
L-Code L4396/L4397
Addresses plantar flexion, flaccid foot drop, plantar fasciitis and mild abnormal muscle tone.



RM Safe Boot I
Product # 20744I
L-Code L4398
Addresses flaccid foot drop and suspends the heel.



RM Fluidized Boot
Product # 20743
L4396/4397
Helps to maintain neutral alignment of the foot and to relieve pressure on the heel using fluidized technology.



Complete Cair Boot
Product # RM20747U
L-Code L4396/L4397
Helps maintain neutral alignment of the foot and to relieve pressure on the heel.



NeuroFlex® Restorative Flex Boot
Product # 20734
L-Code L4396/4397
Addresses mild to severe lost range and abnormal muscle tone.



EZ Stride Anterior Carbon Fiber AFO
Product # CFA
L-Code L1932
Enhances normal gait for persons with flaccid foot drop. The Lateral strut design spirals to anterior aspect, providing support for patients who pronate (valgus or eversion). Tibial pad available.



EZ Stride Posterior Carbon Fiber AFO
Product # CFP
L-Code L1951
Enhances normal gait for persons with flaccid foot drop. The lateral strut design spirals to medial aspect. Best for patients who cannot tolerate tibial pressure.

ELBOW/SHOULDER



NeuroFlex® Restorative Elbow
Product # 20724
L-Code L3760
Addresses mild to severe flexion contractures. Available in styles to treat extensor tone/hyperextended elbows, and for below elbow amputations (BEA).



Elbow Separator (Accessory for NeuroFlex® Restorative Elbow)
Product # 20724SP
Helps to gradually separate the elbow joint before transitioning to the NeuroFlex® Restorative Elbow.



RM RestAir Elbow
Product # 20724air
Suggested L-Code L3760
Addresses elbow flexion contractures with air technology.



NeuroFlex® Restorative Hyper-Extension Elbow
Product # 20724H
L-Code L3760
Helps to improve and/or restore extension lost range of the elbow



Restorative Shoulder
Product # RM6001
L-Code L3960
Provides stability and support for chronic shoulder pain.

RESTORATIVE MEDICAL PRODUCT SUMMARY

Our comprehensive selection of static splints incorporate dynamic properties and technologies to provide therapeutic positioning that will aid in the management and rehabilitation of contractures. Sizes range from infant to adult.

BACK/NECK



RM VertebrEase TLSO
Product # RM20730
L-Code L0456/L0457
Helps improve posture for spinal deformities such as: kyphosis, lordosis, scoliosis and more. Padded shoulder straps encourage upright posture.



RM VertebrEase LSO
Product # 20729
L-Code L0631/L0648
Helps immobilize and support the lumbar spine for soft tissue injuries and relieve lower back pain. Shown with optional shoulder straps. Available with pendulous abdomen front.



NeuroFlex® Restorative Safe Spine
Product # 20728
L-Code L0454/L0455
Comfortable spine support for conditions like kyphosis, scoliosis, lordosis. Padded shoulder straps encourage upright posture.



NeuroFlex® Restorative Kentucky Kollar
Product # 20735
L-Code L0113
Addresses dystonia, torticollis and wry neck conditions.



Thermal Back LSO
Product # RM502014
Lcode L0637/L0650
Reduce inflammation and pain with **cold/warm** pack and adjustable compression.



Restorative LSO 637
Product # RM20637
L-Code L0637/L0650
Helps to improve posture and lumbar spinal stability for herniated disc, stenosis, spondylosis, degenerative disc disease, facet joint syndrome, and post-operative rehabilitation.



Restorative TLSO 457
Product # RM20456
L-Code L0457
Helps to improve posture and spinal stability for herniated disc, stenosis, spondylosis, degenerative disc disease, facet joint syndrome, and post-operative rehabilitation.



Restorative Lumbar Sacral Orthosis (LSO) Stringback
Product # RM10037
Provides support, helps manage trunk stabilization, and improves posture and torso alignment



NeuroFlex® Restorative Cervical HyperExtension Stop
Product # 20728H
L-Code L0999
Fits into our TLSO and addresses neck hyper-extension, head leaning.



NeuroFlex® Restorative Safe Cerv, Left or Right
Product # 20728SCL/R
L-Code L0999
Fits into our TLSO and addresses left or right head leaning with or without forward leaning.



NeuroFlex® Restorative Cervical Extension, Forward
Product # 20728C
L-Code L0999
Fits into our VertebrEase and Safe Spine TLSOs and addresses forward head upper body posture.

KNEE/HIP



NeuroFlex® Restorative Knee
Product # 20725
L-Code L1831
Addresses mild to severe flexion contractures. Available in styles to treat extensor tone/hyperextended knees, and for below knee amputations (BKA).



Knee Separator (Accessory to NeuroFlex® Restorative Knee)
Product # 20725SP
Helps to gradually separate the knee joint before transitioning to the NeuroFlex® Restorative Knee.



RM RestAir Knee
Product # 20755air
L-Code L1831/L1847
Addresses knee flexion contractures with air technology.



NeuroFlex® Restorative Hyperextension Knee
Product # 20725H
Suggested L-Code L1831
Helps to improve and/or restore extension lost range of the knee with effective prolonged low load passive stretch.



RM RestAir Hip
Product # 20760air
L-Code L1652
Comfortably and safely helps to address hip adduction & scissoring.



Restorative OA Knee
Product # 9021
L-Code: L1843
Provides stability and support for the knee.



Restorative Hinged Knee
Product # 9001
L-Code: L1832/L1833
Provides stability and support for chronic knee pains and instability.