

2017

PRODUCT GUIDE

BACK / NECK		NeuroFlex® Restorative Kentucky Kollar	RM VertebrEase TLSO	NeuroFlex® Restorative Safe Spine TLSO	RM VertebrEase LSO	NeuroFlex® Restorative Hyperextension Stop (for VertebrEase or Safe Spine TLSO)	NeuroFlex® Restorative Cervical Extension, forward (for VertebrEase or Safe Spine TLSO)	NeuroFlex® Restorative Safe Cerv, Left or Right (for VertebrEase or Safe Spine TLSO)
								
Item #		#20735	#20730	#20728	#20729	#20728H	#20728C	#20728SC(L or R)
HEAD / NECK	Wry Neck / Torticollis							
	Hyper-extended Head Leaning							
	Left or Right Head Leaning							
	Upper Body/Head Forward							
	Dystonic Conditions with Forward Head Leaning							
Lower Back Pain								
Torso								
Motor or Movement Disorders Related To:								
Multiple Sclerosis, Dementia, ALS, Parkinson's								
TBI, CVA, Stroke, Whiplash								
Cerebral Palsy, End Stage Alzheimer's								

* This chart is intended to provide caregivers a snapshot of different product features. It may not be comprehensive and should not be used as the sole factor in product selection and use. Each Individual's need are unique. A licensed healthcare practitioner should always verify that the product is appropriate for the individual prior to use.

NeuroFlex™ Restorative Kentucky Kollar

Measuring Guide


LCODE: L0113

Item Code	Size	Forehead to Shoulder
20735XXXS	3X Small	2" to 3"
20735XXS	2X Small	2 1/2" to 3 1/2"
20735XS	X Small	3" to 4"
20735S	Small	4" to 5"
20735M	Medium	5" to 6"
20735L	Large	6" to 7"
20735XL	X Large	7" to 9"



2017

PRODUCT GUIDE



*Choose the greater of the two circumference measurements when using the measuring guide.

*For all available sizes please consult our website: www.restorativemedical.com

RM VertebrEase TLSO

Measuring Guide

LCODE: L0456/L0457

Item Code	Size	Back Panel	Front Panel	Waist/Hip Circ.
20730I	Infant	7"	5 1/2"	15-20"
20730T	Toddler	10"	8"	20-25"
20730C	Child	12"	10"	23-30"
20730XSE	X-Short X-Slim	12"	6-7 1/2"	25-30"
20730XSS	X-Short Slim	12"	6-7 1/2"	30-35"
20730XSR	X-Short Regular	12"	6-7 1/2"	35-40"
20730XSF	X-Short Full	12"	6-7 1/2"	40-46"
20730XSX	X-Short X-Full	12"	6-7 1/2"	46-53"
20730SE	Short X-Slim	14"	7-8 1/2"	25-30"
20730SS	Short Slim	14"	7-8 1/2"	30-35"
20730SR	Short Regular	14"	7-8 1/2"	35-40"
20730SF	Short Full	14"	7-8 1/2"	40-46"
20730SX	Short X-Full	14"	7-8 1/2"	46-53"
20730AE	Average X-Slim	16"	8-9 1/2"	25-30"
20730AS	Average Slim	16"	8-9 1/2"	30-35"
20730AR	Average Regular	16"	8-9 1/2"	35-40"
20730AF	Average Full	16"	8-9 1/2"	40-46"
20730AX	Average X-Full	16"	8-9 1/2"	46-53"
20730TE	Tall X-Slim	18"	9-10 1/2"	25-30"
20730TS	Tall Slim	18"	9-10 1/2"	30-35"
20730TR	Tall Regular	18"	9-10 1/2"	35-40"
20730TF	Tall Full	18"	9-10 1/2"	40-46"
20730TX	Tall X-Full	18"	9-10 1/2"	46-53"



NeuroFlex™ Restorative Safe Spine

Measuring Guide

LCODE: L0454/L0455

Item Code	Size	Back Panel	Front Panel	Waist/Hip Circ.
20728I	Infant	7"	5 1/2"	15-20"
20728T	Toddler	10"	8"	20-25"
20728C	Child	12"	10"	23-30"
20728XSE	X-Short X-Slim	12"	6-7 1/2"	25-30"
20728XSS	X-Short Slim	12"	6-7 1/2"	30-35"
20728XSR	X-Short Regular	12"	6-7 1/2"	35-40"
20728XSF	X-Short Full	12"	6-7 1/2"	40-46"
20728XSX	X-Short X-Full	12"	6-7 1/2"	46-53"
20728SE	Short X-Slim	14"	7-8 1/2"	25-30"
20728SS	Short Slim	14"	7-8 1/2"	30-35"
20728SR	Short Regular	14"	7-8 1/2"	35-40"
20728SF	Short Full	14"	7-8 1/2"	40-46"
20728SX	Short X-Full	14"	7-8 1/2"	46-53"
20728AE	Average X-Slim	16"	8-9 1/2"	25-30"
20728AS	Average Slim	16"	8-9 1/2"	30-35"
20728AR	Average Regular	16"	8-9 1/2"	35-40"
20728AF	Average Full	16"	8-9 1/2"	40-46"
20728AX	Average X-Full	16"	8-9 1/2"	46-53"
20728TE	Tall X-Slim	18"	9-10 1/2"	25-30"
20728TS	Tall Slim	18"	9-10 1/2"	30-35"
20728TR	Tall Regular	18"	9-10 1/2"	35-40"
20728TF	Tall Full	18"	9-10 1/2"	40-46"
20728TX	Tall X-Full	18"	9-10 1/2"	46-53"



2017

PRODUCT GUIDE

RM VertebrEase LSO

Measuring Guide

LCODE: **L0631/L0648**

Item Code	Size	Back Panel	Front Panel	Waist/Hip Circ.
20729XSE	X-Short X-Slim	11"	6-7 1/2"	25-30"
20729XSS	X-Short Slim	11"	6-7 1/2"	30-35"
20729XSR	X-Short Regular	11"	6-7 1/2"	35-40"
20729XSF	X-Short Full	11"	6-7 1/2"	40-46"
20729XSX	X-Short X-Full	11"	6-7 1/2"	46-53"
20729SE	Short X-Slim	12"	7-8 1/2"	25-30"
20729SS	Short Slim	12"	7-8 1/2"	30-35"
20729SR	Short Regular	12"	7-8 1/2"	35-40"
20729SF	Short Full	12"	7-8 1/2"	40-46"
20729SXF	Short X-Full	12"	7-8 1/2"	46-53"
20729AE	Average X-Slim	13"	8-9 1/2"	25-30"
20729AS	Average Slim	13"	8-9 1/2"	30-35"
20729AR	Average Regular	13"	8-9 1/2"	35-40"
20729AF	Average Full	13"	8-9 1/2"	40-46"
20729AX	Average X-Full	13"	8-9 1/2"	46-53"
20729TE	Tall X-Slim	14"	9-10 1/2"	25-30"
20729TS	Tall Slim	14"	9-10 1/2"	30-35"
20729TR	Tall Regular	14"	9-10 1/2"	35-40"
20729TF	Tall Full	14"	9-10 1/2"	40-46"
20729TX	Tall X-Full	14"	9-10 1/2"	46-53"



NeuroFlex® Restorative Hyperextension Stop

Measuring Guide

Item Code	Size
20728H	Universal

(for VertebrEase or Safe Spine TLSO)

LCODE: **L0999**



NeuroFlex® Restorative Cervical Extension, Forward

Measuring Guide

Item Code	Size
20728C	Universal

(for VertebrEase or Safe Spine TLSO)

LCODE: **L0999**



NeuroFlex® Restorative Safe Cerv, Left or Right

Measuring Guide

Item Code	Size	Head Circumference
20728SCL-L	Left Large	22" or Greater
20728SCR-L	Right Large	22" or Greater
20728SCL-S	Left Small	Less Than 22"
20728SCR-S	Right Small	Less Than 22"

(for VertebrEase or Safe Spine TLSO)

LCODE: **L0999**




2017

PRODUCT GUIDE

HAND		RM Restorative Hand	RM Prosperity Hand	NeuroFlex® Restorative Flex Hand	RM Thumb Spica	NeuroFlex® Restorative ThumbEase	RM Cock-Up Wrist Splint	RM Cock-Up Hand Splint with Finger Extension Roll	RM BendEase Hand	RM Pediatric HyperHand
Item #		#20736	#20732	#20721	#20719	#20737 / #20740	#20739	#20739cone	#20715	#20717
Hand / Wrist / Forearm	Mild/Moderate Contracture									
	Severe Contracture			Moderate to Closed Palm		Severe Thumb				
Pediatric Sizes Available										
Finger	Severe Flexion in 4th & 5th Digit									
	Severe Thumb ADD									
Moldable By Hand										
Motor or Movement Disorders Related To:										
Cerebral Palsy										
Multiple Sclerosis										
Arthritis										
CVA or Stroke										
Carpal Tunnel Syndrome		Remold Wrist	Remold Wrist	Remold Wrist		Remold Wrist	Molded to 30 degrees	Molded to 30 degrees		

* This chart is intended to provide caregivers a snapshot of different product features. It may not be comprehensive and should not be used as the sole factor in product selection and use. Each Individual's need are unique. A licensed healthcare practitioner should always verify that the product is appropriate for the individual prior to use.



RM Restorative Hand

Measuring Guide

LCODE: **L3807/3809**

Item Code	Size	Width of MCP	Wrist to Finger
20736XSR	Ex Small -Right	2 3/4" to 3"	Up to 5"
20736XSL	Ex Small -Left	2 3/4" to 3"	Up to 5"
20736SR	Small -Right	Up to 3"	Up to 6 1/2"
20736SL	Small -Left	Up to 3"	Up to 6 1/2"
20736MR	Medium -Right	Up to 3 1/2"	Up to 7"
20736ML	Medium -Left	Up to 3 1/2"	Up to 7"
20736LR	Large -Right	Up to 3 3/4"	Up to 8"
20736LL	Large -Left	Up to 3 3/4"	Up to 8"
20736RZ	Custom - Right		
20736LZ	Custom - Left		



RM Prosperity Hand

Measuring Guide

LCODE: L3807/3809

Item Code	Size	Width of MCP	Wrist to Finger
20732SR	Small -Right	Up to 3"	Up to 6 1/2"
20732SL	Small -Left	Up to 3"	Up to 6 1/2"
20732MR	Medium -Right	Up to 3 1/2"	Up to 7"
20732ML	Medium -Left	Up to 3 1/2"	Up to 7"
20732LR	Large -Right	Up to 3 3/4"	Up to 8"
20732LL	Large -Left	Up to 3 3/4"	Up to 8"



NeuroFlex® Restorative Flex Hand

Measuring Guide

LCODE: L3807/3809

Item Code	Size	Width of MCP
20720TR	Toddler -Right	2" to 2 1/4"
20720TL	Toddler -Left	2" to 2 1/4"
20720CR	Child-Right	2 1/4" to 2 1/2"
20720CL	Child-Left	2 1/4" to 2 1/2"
20720YR	Youth -Right	2 3/4" to 3"
20720YL	Youth -Left	2 3/4" to 3"

Universal Adult Size LF or RT with cones or cylinders)

20721L	Left w/Regular Cones
20721R	Right w/Regular Cones
20721LRC	Left w/Reverse Cones
20721RRC	Right w/Reverse Cones
20721LCY	Left w/Cylinders
20721RCY	Right w/Cylinders
20721LZ	Custom C-BAR Left
20721RZ	Custom C-BAR Right




RM Thumb Spica

Measuring Guide

LCODE: L3807/3809

Item Code	Size	Width of MCP
20719IR	Infant -Right	1 3/4" to 2"
20719IL	Infant -Left	1 3/4" to 2"
20719TJR	Toddler Jr. -Right	1 3/4" to 2"
20719TJL	Toddler Jr. -Left	1 3/4" to 2"
20719TR	Toddler -Right	2" to 2 1/4"
20719TL	Toddler -Left	2" to 2 1/4"
20719CR	Child-Right	2 1/4" to 2 1/2"
20719CL	Child-Left	2 1/4" to 2 1/2"
20719XSR	Ex Small -Right	2 3/4" to 3"
20719XSL	Ex Small -Left	2 3/4" to 3"
20719SR	Small -Right	Up to 3"
20719SL	Small -Left	Up to 3"
20719MR	Medium -Right	Up to 3 1/2"
20719ML	Medium -Left	Up to 3 1/2"
20719LR	Large -Right	Up to 3 3/4"
20719LL	Large -Left	Up to 3 3/4"



NeuroFlex® Restorative ThumbEase (Pediatric)

Measuring Guide

LCODE: L3807/3809

Item Code	Size	Width of MCP	Wrist to Finger	Wrist to Forearm
20740TR	Toddler -Right	2" to 2 1/4"	Up To 3 1/2"	3 1/2"
20740TL	Toddler -Left	2" to 2 1/4"	Up To 3 1/2"	3 1/2"
20740CR	Child -Right	2 1/4" to 2 1/2"	Up To 4"	3 3/4"
20740CL	Child -Left	2 1/4" to 2 1/2"	Up To 4"	3 3/4"
20740YR	Youth -Right	2 3/4" to 3"	Up To 5"	4"
20740YL	Youth -Left	2 3/4" to 3"	Up To 5"	4"



NeuroFlex® Restorative ThumbEase (Adult)

Measuring Guide

LCODE: L3807/3809

Item Code	Size	Width of MCP	Wrist to Finger
20737XSR	Ex Small -Right	2 3/4" to 3"	Up to 5"
20737XSL	Ex Small -Left	2 3/4" to 3"	Up to 5"
20737SR	Small -Right	Up to 3"	Up to 6 1/2"
20737SL	Small -Left	Up to 3"	Up to 6 1/2"
20737MR	Medium -Right	Up to 3 1/2"	Up to 7"
20737ML	Medium -Left	Up to 3 1/2"	Up to 7"
20737LR	Large -Right	Up to 3 3/4"	Up to 8"
20737LL	Large -Left	Up to 3 3/4"	Up to 8"
20737RZ	Custom -Right		
20737LZ	Custom -Left		



RM Cock-Up Wrist Splint

Measuring Guide

Item Code	Size	Width of MCP
20739XSR	Ex Small -Right	2 3/4" to 3"
20739XSL	Ex Small -Left	2 3/4" to 3"
20739SR	Small -Right	Up to 3"
20739SL	Small -Left	Up to 3"
20739MR	Medium -Right	Up to 3 1/2"
20739ML	Medium -Left	Up to 3 1/2"
20739LR	Large -Right	Up to 3 3/4"
20739LL	Large -Left	Up to 3 3/4"



RM Cock-Up Hand Splint with Finger Extension Roll

Measuring Guide

LCODE: L3807/3809

Item Code	Size	Width of MCP
20739rollXSR	Ex Small -Right	2 3/4" to 3"
20739rollXSL	Ex Small -Left	2 3/4" to 3"
20739rollISR	Small -Right	Up to 3"
20739rollISL	Small -Left	Up to 3"
20739rollIMR	Medium -Right	Up to 3 1/2"
20739rollIML	Medium -Left	Up to 3 1/2"
20739rollILR	Large -Right	Up to 3 3/4"
20739rollILL	Large -Left	Up to 3 3/4"



RM BendEase Hand

Measuring Guide

LCODE: L3807/3809

Item Code	Size	Width of MCP	Wrist to Finger
20715SR	Small -Right	Up to 3"	Up to 6 1/2"
20715SL	Small -Left	Up to 3"	Up to 6 1/2"
20715MR	Medium -Right	Up to 3 1/2"	Up to 7"
20715ML	Medium -Left	Up to 3 1/2"	Up to 7"
20715LR	Large -Right	Up to 3 3/4"	Up to 8"
20715LL	Large -Left	Up to 3 3/4"	Up to 8"



RM BendEase Hand (Universal, Left or Right)

Measuring Guide

LCODE: L3807/3809

Item Code	Size
20715L	Universal Left
20715R	Universal Right



RM Pediatric HyperHand

Measuring Guide

LCODE: L3807/3809

Item Code	Size	Width of MCP	Wrist to Finger	Wrist to Forearm
20717IR	Infant -Right	1 3/4" to 2"	Up To 2"	2 1/4"
20717IL	Infant -Left	1 3/4" to 2"	Up To 2"	2 1/4"
20717TJR	Toddler Jr. -Right	1 3/4" to 2"	Up To 3"	2 3/4"
20717TJL	Toddler Jr. -Left	1 3/4" to 2"	Up To 3"	2 3/4"
20717TR	Toddler -Right	2" to 2 1/4"	Up To 3 1/2"	3 1/2"
20717TL	Toddler -Left	2" to 2 1/4"	Up To 3 1/2"	3 1/2"
20717CR	Child -Right	2 1/4" to 2 1/2"	Up To 4"	3 3/4"
20717CL	Child -Left	2 1/4" to 2 1/2"	Up To 4"	3 3/4"
20717YR	Youth -Right	2 3/4" to 3"	Up To 5"	4"
20717YL	Youth -Left	2 3/4" to 3"	Up To 5"	4"



Accessories - Hand Splint

20736D	Dorsal Interphalangeal Flexion Assist
20736IS	Infinity Dorsal Attachment - Small
20736IM	Infinity Dorsal Attachment - Medium
20736IL	Infinity Dorsal Attachment - Large
20736V	Ventral Interphalangeal Flexion Assist "ROLL"
20737FS	NeuroFlex® Restorative™ Velcro-on Finger Separators
20737UD	NeuroFlex® Restorative™ Velcro-on Ulnar Drift Finger Separators
20737UDstraps	Sewn-on RMI Ulnar Drift Wrist Straps - addition available on hand splint



20737UD -
Ulnar Drift Finger Separators



20737UDstraps -
Ulnar Drift Strapping



20736D -
Dorsal Interphalangeal Flexion Assist



20736I -
Infinity Dorsal Attachment

ACCESSORIES


2017

PRODUCT GUIDE

<i>ELBOW / KNEE / HIP</i>		NeuroFlex® Restorative Elbow	NeuroFlex® Restorative HyperExtension Elbow	RM RestAir Elbow	NeuroFlex® Restorative Knee	NeuroFlex® Restorative HyperExtension Knee	RM RestAir Knee	RM RestAir Hip
Item #		#20724	#20724H	#20724air	#20725	#20725H	#20755	#20760
Elbow	Mild/Moderate Contracture							
	Severe Contracture							
Knee	Mild/Moderate Contracture							
	Severe Contracture							
Hip								Adduction/Scissoring
Motor or Movement Disorders Related To:								
CVA or Stroke								
Rheumatoid Arthritis								
Osteoarthritis								
TBI, ABI, Stroke								
Parkinson's, Dementia, End Stage Alzheimer's								
Extensor Tone/Hyperextension								
Cerebral Palsy, Multiple Sclerosis								
* This chart is intended to provide caregivers a snapshot of different product features. It may not be comprehensive and should not be used as the sole factor in product selection and use. Each individual's needs are unique. A licensed healthcare practitioner should always verify that the product is appropriate for the individual prior to use.								

2017

PRODUCT GUIDE



NeuroFlex® Restorative Elbow

Measuring Guide

LCODE: L3760

Item Code	Size	BTE	ETF	Biceps Cir.	Forearm Cir.
20724C	Child	3 ½"	4"	5 ½ - 7 ½"	3 ½ - 6 ½"
20724CF	Child/Full	3 ½"	4"	6 ½ - 8 ½"	5 ½ - 7 ½"
20724YS	Youth/Slim	4"	5"	5 ½ - 7 ½"	3 ½ - 6 ½"
20724Y	Youth	4"	5"	6 ½ - 8 ½"	5 ½ - 7 ½"
20724XS	Ex-Small	5"	5"	7 ½ - 9 ½"	6 ½ - 8 ½"
20724XSF	Ex-Small/Full	5"	5"	8-12"	7-10"
20724SS	Small/Slim	5"	6"	7 ½ - 9 ½"	6 ½ - 8 ½"
20724S	Small	5"	6"	8-12"	7-10"
20724SF	Small/Full	5"	6"	8-13"	7 ½ - 11"
20724MS	Medium/Slim	6"	6"	8-12"	7-10"
20724M	Medium	6"	6"	8-13"	7 ½ - 11"
20724MF	Medium/Full	6"	6"	10-15"	8-12"
20724LS	Large/Slim	6 ½"	6 ½"	8-13"	7 ½ - 11"
20724L	Large	6 ½"	6 ½"	10-15"	8-12"
20724LF	Large/Full	6 ½"	6 ½"	13-17"	10-15"
20724Z	Custom				



NeuroFlex® Restorative HyperExtension Elbow

Measuring Guide

LCODE: L3760

Item Code	Size	BTE	ETF	Biceps Cir.	Forearm Cir.
20724CH	Child	3 ½"	4"	5 ½ - 7 ½"	3 ½ - 6 ½"
20724CFH	Child/Full	3 ½"	4"	6 ½ - 8 ½"	5 ½ - 7 ½"
20724YSH	Youth/Slim	4"	5"	5 ½ - 7 ½"	3 ½ - 6 ½"
20724YH	Youth	4"	5"	6 ½ - 8 ½"	5 ½ - 7 ½"
20724XSH	Ex-Small	5"	5"	7 ½ - 9 ½"	6 ½ - 8 ½"
20724XSFH	Ex-Small/Full	5"	5"	8-12"	7-10"
20724SSH	Small/Slim	5"	6"	7 ½ - 9 ½"	6 ½ - 8 ½"
20724SH	Small	5"	6"	8-12"	7-10"
20724SFH	Small/Full	5"	6"	8-13"	7 ½ - 11"
20724MSH	Medium/Slim	6"	6"	8-12"	7-10"
20724MH	Medium	6"	6"	8-13"	7 ½ - 11"
20724MFH	Medium/Full	6"	6"	10-15"	8-12"
20724LSH	Large/Slim	6 ½"	6 ½"	8-13"	7 ½ - 11"
20724LH	Large	6 ½"	6 ½"	10-15"	8-12"
20724LFH	Large/Full	6 ½"	6 ½"	13-17"	10-15"
20724ZH	Custom				



RM RestAir Elbow

Measuring Guide

LCODE: L3760

Item Code	Arm	Total Length	Upper Arm Circ.	Lower Arm Circ.
20724airR	Right	12"	8-15"	7-12"
20724airL	Left	12"	8-15"	7-12"



NeuroFlex® Restorative Elbow Separator


Measuring Guide

Item Code	Size
20724SP	Universal-Adult



2017

PRODUCT GUIDE



NeuroFlex® Restorative Knee

Measuring Guide

LCODE: L1831

Item Code	Size	TTK	KTC	Thigh Cir.	Calf Cir.
20725CS	Child/Slim	5"	6"	7 ½ - 9 ½ "	6 ½ - 8 ½ "
20725C	Child	5"	6"	8-12"	7-10"
20725CF	Child/Full	5"	6"	10-15"	8-12"
20725XSS	Ex-Small/Slim	6"	6"	8-12"	7-10"
20725XS	Ex-Small*	6"	6"	8-13"	7 1/2-11"
20725XSF	Ex-Small/Full	6"	6"	10-15"	8-12"
20725SS	Small/Slim	6 ½ "	6 ½ "	8-13"	7 1/2-11"
20725S	Small	6 ½ "	6 ½ "	10-15"	8-12"
20725SF	Small/Full	6 ½ "	6 ½ "	13-17"	10-15"
20725MS	Medium/Slim	7 ½ "	7 ½ "	10-15"	8-12"
20725M	Medium	7 ½ "	7 ½ "	13-17"	10-15"
20725MF	Medium/Full	7 ½ "	7 ½ "	16-22"	14-18"
20725LS	Large/Slim	8 ½ "	8 ½ "	13-17"	10-15"
20725L	Large	8 ½ "	8 ½ "	16-22"	14-18"
20725Z	Custom				



NeuroFlex® Restorative HyperExtension Knee

Measuring Guide

LCODE: L1831

Item Code	Size	TTK	KTC	Thigh Cir.	Calf Cir.
20725CSH	Child/Slim	5"	6"	7 ½ - 9 ½ "	6 ½ - 8 ½ "
20725CH	Child	5"	6"	8-12"	7-10"
20725CFH	Child/Full	5"	6"	10-15"	8-12"
20725XSSH	Ex-Small/Slim	6"	6"	8-12"	7-10"
20725XSH	Ex-Small	6"	6"	8-13"	7 1/2-11"
20725XSFH	Ex-Small/Full	6"	6"	10-15"	8-12"
20725SSH	Small/Slim	6 ½ "	6 ½ "	8-13"	7 1/2-11"
20725SH	Small	6 ½ "	6 ½ "	10-15"	8-12"
20725SFH	Small/Full	6 ½ "	6 ½ "	13-17"	10-15"
20725MSH	Medium/Slim	7 ½ "	7 ½ "	10-15"	8-12"
20725MH	Medium	7 ½ "	7 ½ "	13-17"	10-15"
20725MFH	Medium/Full	7 ½ "	7 ½ "	16-22"	14-18"
20725LSH	Large/Slim	8 ½ "	8 ½ "	13-17"	10-15"
20725LH	Large	8 ½ "	8 ½ "	16-22"	14-18"
20725ZH	Custom				



NeuroFlex® Restorative Knee Separator

Measuring Guide

Item Code	Size
20725SP	Universal-Adult



RM RestAir Knee

Measuring Guide

LCODE: L1847

Item Code	Size	Total Length	Thigh Circ.	Calf Circ.
20755SR	Small -Right	15"	10-15"	8-12"
20755SL	Small-Left	15"	10-15"	8-12"
20755MR	Medium-Right	15"	13-17"	10-15"
20755ML	Medium-Left	15"	13-17"	10-15"
20755LR	Large-Right	17"	17-22"	12-17"
20755LL	Large-Left	17"	17-22"	12-17"



RM RestAir Hip

Measuring Guide

LCODE: L1652

Item Code	Size	Thigh Circ
20760airXS	X-Small	6"-12"
20760airS	Small	11"-17"
20760airL	Large	16"-22"
20760airXL	X-Large	21"-26"




2017

PRODUCT GUIDE

<i>Foot / Ankle</i>		RM Pediatric Safe Boot	NeuroFlex® Restorative Flex Boot	NeuroFlex® Restorative FootEase Boot	RM Safe Boot II	RM Safe Boot I	RM Flo-Form Ankle Foot Orthosis	EZ Stride Anterior Carbon Fiber AFO	EZ Stride Posterior Carbon Fiber AFO
								 *Distributed by Restorative Medical	 *Distributed by Restorative Medical
Item #		#20734	#20734	#20742	#20744I	#20744II	SUN TAFO	CFA	CFP
Foot / Ankle	Mild/Moderate Contracture								
	Severe Contracture								
	Foot Drop					Flaccid Foot Drop		Flaccid Foot Drop	Flaccid Foot Drop
	Heel Suspension / Protection								
	Plantar Fasciitis								
Motor or Movement Disorders Related To:									
CVA, Stroke, TBI, ABI									
Rheumatoid Arthritis									
Osteoarthritis									
ALS									
Spinal Cord Injuries									
Parkinson's, Dementia, End Stage Alzheimer's, Cerebral Palsy									

* This chart is intended to provide caregivers a snapshot of different product features. It may not be comprehensive and should not be used as the sole factor in product selection and use. Each individual's need are unique. A licensed healthcare practitioner should always verify that the product is appropriate for the individual prior to use.



2017

PRODUCT GUIDE

RM Pediatric Safe Boot

Measuring Guide

LCODE: L4396/L4397

Item Code	Size	Foot Length	Infant Shoe	Kids Shoe	Youth Shoe	Men's Shoe	Women's Shoe
20746-3	Pediatric 3"	3"	0.5				
20746-3.5	Pediatric 3.5"	3.5"	1				
20746-4	Pediatric 4"	4"	2.5				
20746-5	Pediatric 5"	5"		5.5			
20746-6	Pediatric 6"	6"		8.5			
20746-7	Pediatric 7"	7"		11.5			
20746-8	Pediatric 8"	8"		13			
20746Z	Custom						



NeuroFlex® Restorative Flex Boot

Measuring Guide

LCODE: L4396/L4397

Item Code	Size	Foot Length	Infant Shoe	Kids Shoe	Youth Shoe	Men's Shoe	Women's Shoe
20734-S	Small	Up to 8"			1.5		4
20734-M	Medium	Up to 9 1/2 "				6.5	8
20734-L	Large	Up to 11"				11	12
20734-XL	Ex-Large	Up to 13"				16	
20734Z	Custom						



NeuroFlex® Restorative FootEase Boot

Measuring Guide

LCODE: L4396/L4397

Item Code	Size	Foot Length	Infant Shoe	Kids Shoe	Youth Shoe	Men's Shoe	Women's Shoe
20742-S	Small	Up to 8"			1.5		4
20742-M	Medium	Up to 9 1/2 "				6.5	8
20742-L	Large	Up to 11"				11	12
20742-XL	Ex-Large	Up to 13"				16	



RM Safe Boot I

Measuring Guide

LCODE: L4398

Item Code	Size	Foot Length	Infant Shoe	Kids Shoe	Youth Shoe	Men's Shoe	Women's Shoe
20744I	Universal						
20744I-S	Small	Up to 8"			1.5		4
20744I-M	Medium	Up to 9 1/2 "				6.5	8
20744I-L	Large	Up to 11"				11	12
20744I-XL	Ex-Large	Up to 13"				16	



RM Safe Boot II

Measuring Guide

LCODE: L4396/L4397

Item Code	Size	Foot Length	Infant Shoe	Kids Shoe	Youth Shoe	Men's Shoe	Women's Shoe
20744II	Universal						
20744II-S	Small	Up to 8"			1.5		4
20744II-M	Medium	Up to 9 1/2 "				6.5	8
20744II-L	Large	Up to 11"				11	12
20744II-XL	Ex-Large	Up to 13"				16	



Fluidized Ankle Foot Orthosis

Measuring Guide

LCODE: L4396/L4397

Item Code	Size	Foot Length	Infant Shoe	Kids Shoe	Youth Shoe	Men's Shoe	Women's Shoe
SUN TAFO -S	Small	Up to 8"			1.5		4
SUN TAFO -M	Medium	Up to 9 1/2 "				6.5	8
SUN TAFO -L	Large	Up to 11"				11	12
SUN TAFO -XL	Ex-Large	Up to 13"				16	



RM Flo-Form Complete Care Boot

Measuring Guide

LCODE: L4396/L4397

Item Code	Size	Foot Length	Infant Shoe	Kids Shoe	Youth Shoe	Men's Shoe	Women's Shoe
20747-S	Small	Up to 8"			1.5		4
20747-M	Medium	Up to 9 1/2 "				6.5	8
20747-L	Large	Up to 11"				11	12
20747-XL	Ex-Large	Up to 13"				16	

EZ Stride Anterior Carbon Fiber AFO

Measuring Guide

LCODE: **L1932**

Item Code	Size	Footplate Length	Height Length
CFA-3XSR	XXX-Small Right	7	13"
CFA-3XSL	XXX-Small Left	7	13"
CFA-2XSR	XX-Small Right	8	14"
CFA-2XSL	XX-Small Left	8	14"
CFA-XSR	Ex- Small Right	9"	14
CFA-XSL	Ex- Small Left	9"	14
CFA-SR	Small Right	10"	15
CFA-SL	Small Left	10"	15
CFA-MR	Medium Right	11	16
CFA-ML	Medium Left	11	16
CFA-LR	Large Right	12	17"
CFA-LL	Large Left	12	17"
CFA-Kit	Universal Pad Kit	<i>Suggested code L2820</i>	

* DISTRIBUTED BY RESTORATIVE MEDICAL



EZ Stride Posterior Carbon Fiber AFO

Measuring Guide

LCODE: **L1951**

Item Code	Size	Footplate Length	Height Length
CFP-3XSR	XXX-Small Right	7	13"
CFP-3XSL	XXX-Small Left	7	13"
CFP-2XSR	XX-Small Right	8	14"
CFP-2XSL	XX-Small Left	8	14"
CFP-XSR	Ex- Small Right	9"	14
CFP-XSL	Ex- Small Left	9"	14
CFP-SR	Small Right	10"	15
CFP-SL	Small Left	10"	15
CFP-MR	Medium Right	11	16
CFP-ML	Medium Left	11	16
CFP-LR	Large Right	12	17"
CFP-LL	Large Left	12	17"

* DISTRIBUTED BY RESTORATIVE MEDICAL

