




2017 PRODUCT GUIDE

BACK / NECK		NeuroFlex® Restorative Kentucky Kollar	RM VertebrEase TLSO	NeuroFlex® Restorative Safe Spine TLSO	RM VertebrEase LSO	NeuroFlex® Restorative Hyperextension Stop (for VertebrEase or Safe Spine TLSO)	NeuroFlex® Restorative Cervical Extension, forward (for VertebrEase or Safe Spine TLSO)	NeuroFlex® Restorative Safe Cerv, Left or Right (for VertebrEase or Safe Spine TLSO)
								
Item #		#20735	#20730	#20728	#20729	#20728H	#20728C	#207285C(L or R)
HEAD / NECK	Wry Neck / Torticollis							
	Hyper-extended Head Leaning							
	Left or Right Head Leaning							
	Upper Body/Head Forward							
	Dystonic Conditions with Forward Head Leaning							
Lower Back Pain								
Torso								
Motor or Movement Disorders Related To:								
Multiple Sclerosis, Dementia, ALS, Parkinson's								
TBI, CVA, Stroke, Whiplash								
Cerebral Palsy, End Stage Alzheimer's								

* This chart is intended to provide caregivers a snapshot of different product features. It may not be comprehensive and should not be used as the sole factor in product selection and use. Each Individual's need are unique. A licensed healthcare practitioner should always verify that the product is appropriate for the individual prior to use.


NeuroFlex™ Restorative Kentucky Kollar

Measuring Guide

LCODE: L0113

Item Code	Size	Forehead to Shoulder
20735XXXS	3X Small	2" to 3"
20735XXS	2X Small	2 1/2" to 3 1/2"
20735XS	X Small	3" to 4"
20735S	Small	4" to 5"
20735M	Medium	5" to 6"
20735L	Large	6" to 7"
20735XL	X Large	7" to 9"





*Choose the greater of the two circumference measurements when using the measuring guide.

*For all available sizes please consult our website: www.restorativemedical.com

RM VertebrEase TLSO

Measuring Guide

LCODE: **L0456/L0457**

Item Code	Size	Back Panel	Front Panel	Waist/Hip Circ.
20730I	Infant	7"	5 1/2"	15-20"
20730T	Toddler	10"	8"	20-25"
20730C	Child	12"	10"	23-30"
20730XSE	X-Short X-Slim	12"	6-7 1/2"	25-30"
20730XSS	X-Short Slim	12"	6-7 1/2"	30-35"
20730XSR	X-Short Regular	12"	6-7 1/2"	35-40"
20730XSF	X-Short Full	12"	6-7 1/2"	40-46"
20730XSX	X-Short X-Full	12"	6-7 1/2"	46-53"
20730SE	Short X-Slim	14"	7-8 1/2"	25-30"
20730SS	Short Slim	14"	7-8 1/2"	30-35"
20730SR	Short Regular	14"	7-8 1/2"	35-40"
20730SF	Short Full	14"	7-8 1/2"	40-46"
20730SX	Short X-Full	14"	7-8 1/2"	46-53"
20730AE	Average X-Slim	16"	8-9 1/2"	25-30"
20730AS	Average Slim	16"	8-9 1/2"	30-35"
20730AR	Average Regular	16"	8-9 1/2"	35-40"
20730AF	Average Full	16"	8-9 1/2"	40-46"
20730AX	Average X-Full	16"	8-9 1/2"	46-53"
20730TE	Tall X-Slim	18"	9-10 1/2"	25-30"
20730TS	Tall Slim	18"	9-10 1/2"	30-35"
20730TR	Tall Regular	18"	9-10 1/2"	35-40"
20730TF	Tall Full	18"	9-10 1/2"	40-46"
20730TX	Tall X-Full	18"	9-10 1/2"	46-53"



NeuroFlex™ Restorative Safe Spine

Measuring Guide

LCODE: **L0454/L0455**

Item Code	Size	Back Panel	Front Panel	Waist/Hip Circ.
20728I	Infant	7"	5 1/2"	15-20"
20728T	Toddler	10"	8"	20-25"
20728C	Child	12"	10"	23-30"
20728XSE	X-Short X-Slim	12"	6-7 1/2"	25-30"
20728XSS	X-Short Slim	12"	6-7 1/2"	30-35"
20728XSR	X-Short Regular	12"	6-7 1/2"	35-40"
20728XSF	X-Short Full	12"	6-7 1/2"	40-46"
20728XSX	X-Short X-Full	12"	6-7 1/2"	46-53"
20728SE	Short X-Slim	14"	7-8 1/2"	25-30"
20728SS	Short Slim	14"	7-8 1/2"	30-35"
20728SR	Short Regular	14"	7-8 1/2"	35-40"
20728SF	Short Full	14"	7-8 1/2"	40-46"
20728SX	Short X-Full	14"	7-8 1/2"	46-53"
20728AE	Average X-Slim	16"	8-9 1/2"	25-30"
20728AS	Average Slim	16"	8-9 1/2"	30-35"
20728AR	Average Regular	16"	8-9 1/2"	35-40"
20728AF	Average Full	16"	8-9 1/2"	40-46"
20728AX	Average X-Full	16"	8-9 1/2"	46-53"
20728TE	Tall X-Slim	18"	9-10 1/2"	25-30"
20728TS	Tall Slim	18"	9-10 1/2"	30-35"
20728TR	Tall Regular	18"	9-10 1/2"	35-40"
20728TF	Tall Full	18"	9-10 1/2"	40-46"
20728TX	Tall X-Full	18"	9-10 1/2"	46-53"



RM VertebrEase LSO

Measuring Guide

LCODE: **L0631/L0648**

Item Code	Size	Back Panel	Front Panel	Waist/Hip Circ.
20729XSE	X-Short X-Slim	11"	6-7 1/2"	25-30"
20729XSS	X-Short Slim	11"	6-7 1/2"	30-35"
20729XSR	X-Short Regular	11"	6-7 1/2"	35-40"
20729XSF	X-Short Full	11"	6-7 1/2"	40-46"
20729XSX	X-Short X-Full	11"	6-7 1/2"	46-53"
20729SE	Short X-Slim	12"	7-8 1/2"	25-30"
20729SS	Short Slim	12"	7-8 1/2"	30-35"
20729SR	Short Regular	12"	7-8 1/2"	35-40"
20729SF	Short Full	12"	7-8 1/2"	40-46"
20729SXF	Short X-Full	12"	7-8 1/2"	46-53"
20729AE	Average X-Slim	13"	8-9 1/2"	25-30"
20729AS	Average Slim	13"	8-9 1/2"	30-35"
20729AR	Average Regular	13"	8-9 1/2"	35-40"
20729AF	Average Full	13"	8-9 1/2"	40-46"
20729AX	Average X-Full	13"	8-9 1/2"	46-53"
20729TE	Tall X-Slim	14"	9-10 1/2"	25-30"
20729TS	Tall Slim	14"	9-10 1/2"	30-35"
20729TR	Tall Regular	14"	9-10 1/2"	35-40"
20729TF	Tall Full	14"	9-10 1/2"	40-46"
20729TX	Tall X-Full	14"	9-10 1/2"	46-53"



NeuroFlex® Restorative Hyperextension Stop

Measuring Guide

Item Code	Size
20728H	Universal

(for VertebrEase or Safe Spine TLSO)

LCODE: **L0999**



NeuroFlex® Restorative Cervical Extension, Forward

Measuring Guide

Item Code	Size
20728C	Universal

(for VertebrEase or Safe Spine TLSO)

LCODE: **L0999**



NeuroFlex® Restorative Safe Cerv, Left or Right

Measuring Guide

Item Code	Size	Head Circumference
20728SCL-L	Left Large	22" or Greater
20728SCR-L	Right Large	22" or Greater
20728SCL-S	Left Small	Less Than 22"
20728SCR-S	Right Small	Less Than 22"

(for VertebrEase or Safe Spine TLSO)

LCODE: **L0999**

