



### NeuroFlex® Restorative Elbow

#### Measuring Guide

LCODE: L3760

Item Code	Size	BTE	ETF	Biceps Cir.	Forearm Cir.
20724C	Child	3 ½"	4"	5 ½ - 7 ½"	3 ½ - 6 ½"
20724CF	Child/Full	3 ½"	4"	6 ½ - 8 ½"	5 ½ - 7 ½"
20724YS	Youth/Slim	4"	5"	5 ½ - 7 ½"	3 ½ - 6 ½"
20724Y	Youth	4"	5"	6 ½ - 8 ½"	5 ½ - 7 ½"
20724XS	Ex-Small	5"	5"	7 ½ - 9 ½"	6 ½ - 8 ½"
20724XSF	Ex-Small/Full	5"	5"	8-12"	7-10"
20724SS	Small/Slim	5"	6"	7 ½ - 9 ½"	6 ½ - 8 ½"
20724S	Small	5"	6"	8-12"	7-10"
20724SF	Small/Full	5"	6"	8-13"	7 ½ - 11"
20724MS	Medium/Slim	6"	6"	8-12"	7-10"
20724M	Medium	6"	6"	8-13"	7 ½ - 11"
20724MF	Medium/Full	6"	6"	10-15"	8-12"
20724LS	Large/Slim	6 ½"	6 ½"	8-13"	7 ½ - 11"
20724L	Large	6 ½"	6 ½"	10-15"	8-12"
20724LF	Large/Full	6 ½"	6 ½"	13-17"	10-15"
20724Z	Custom				



### NeuroFlex® Restorative HyperExtension Elbow

#### Measuring Guide

LCODE: L3760

Item Code	Size	BTE	ETF	Biceps Cir.	Forearm Cir.
20724CH	Child	3 ½"	4"	5 ½ - 7 ½"	3 ½ - 6 ½"
20724CFH	Child/Full	3 ½"	4"	6 ½ - 8 ½"	5 ½ - 7 ½"
20724YSH	Youth/Slim	4"	5"	5 ½ - 7 ½"	3 ½ - 6 ½"
20724YH	Youth	4"	5"	6 ½ - 8 ½"	5 ½ - 7 ½"
20724XSH	Ex-Small	5"	5"	7 ½ - 9 ½"	6 ½ - 8 ½"
20724XSFH	Ex-Small/Full	5"	5"	8-12"	7-10"
20724SSH	Small/Slim	5"	6"	7 ½ - 9 ½"	6 ½ - 8 ½"
20724SH	Small	5"	6"	8-12"	7-10"
20724SFH	Small/Full	5"	6"	8-13"	7 ½ - 11"
20724MSH	Medium/Slim	6"	6"	8-12"	7-10"
20724MH	Medium	6"	6"	8-13"	7 ½ - 11"
20724MFH	Medium/Full	6"	6"	10-15"	8-12"
20724LSH	Large/Slim	6 ½"	6 ½"	8-13"	7 ½ - 11"
20724LH	Large	6 ½"	6 ½"	10-15"	8-12"
20724LFH	Large/Full	6 ½"	6 ½"	13-17"	10-15"
20724ZH	Custom				



### RM RestAir Elbow

#### Measuring Guide

LCODE: L3760

Item Code	Arm	Total Length	Upper Arm Circ.	Lower Arm Circ.
20724airR	Right	12"	8-15"	7-12"
20724airL	Left	12"	8-15"	7-12"



### NeuroFlex® Restorative Elbow Separator

#### Measuring Guide

Item Code	Size
20724SP	Universal-Adult

