

2017 PRODUCT GUIDE

RM Pediatric Safe Boot

Measuring Guide

LCODE: L4396/L4397

Item Code	Size	Foot Length	Infant Shoe	Kids Shoe	Youth Shoe	Men's Shoe	Women's Shoe
20746-3	Pediatric 3"	3"	0.5				
20746-3.5	Pediatric 3.5"	3.5"	1				
20746-4	Pediatric 4"	4"	2.5				
20746-5	Pediatric 5"	5"		5.5			
20746-6	Pediatric 6"	6"		8.5			
20746-7	Pediatric 7"	7"		11.5			
20746-8	Pediatric 8"	8"		13			
20746Z	Custom						



NeuroFlex® Restorative Flex Boot

Measuring Guide

LCODE: L4396/L4397

Item Code	Size	Foot Length	Infant Shoe	Kids Shoe	Youth Shoe	Men's Shoe	Women's Shoe
20734-S	Small	Up to 8"			1.5		4
20734-M	Medium	Up to 9 1/2"				6.5	8
20734-L	Large	Up to 11"				11	12
20734-XL	Ex-Large	Up to 13"				16	
20734Z	Custom						



NeuroFlex® Restorative FootEase Boot

Measuring Guide

LCODE: L4396/L4397

Item Code	Size	Foot Length	Infant Shoe	Kids Shoe	Youth Shoe	Men's Shoe	Women's Shoe
20742-S	Small	Up to 8"			1.5		4
20742-M	Medium	Up to 9 1/2"				6.5	8
20742-L	Large	Up to 11"				11	12
20742-XL	Ex-Large	Up to 13"				16	



RM Safe Boot I

Measuring Guide

LCODE: L4398

Item Code	Size	Foot Length	Infant Shoe	Kids Shoe	Youth Shoe	Men's Shoe	Women's Shoe
20744I	Universal						
20744I-S	Small	Up to 8"			1.5		4
20744I-M	Medium	Up to 9 1/2"				6.5	8
20744I-L	Large	Up to 11"				11	12
20744I-XL	Ex-Large	Up to 13"				16	



RM Safe Boot II

Measuring Guide

LCODE: L4396/L4397

Item Code	Size	Foot Length	Infant Shoe	Kids Shoe	Youth Shoe	Men's Shoe	Women's Shoe
20744II	Universal						
20744II-S	Small	Up to 8"			1.5		4
20744II-M	Medium	Up to 9 1/2"				6.5	8
20744II-L	Large	Up to 11"				11	12
20744II-XL	Ex-Large	Up to 13"				16	



Fluidized Ankle Foot Orthosis

Measuring Guide

LCODE: L4396/L4397

Item Code	Size	Foot Length	Infant Shoe	Kids Shoe	Youth Shoe	Men's Shoe	Women's Shoe
SUN TAFO -S	Small	Up to 8"			1.5		4
SUN TAFO -M	Medium	Up to 9 1/2"				6.5	8
SUN TAFO -L	Large	Up to 11"				11	12
SUN TAFO -XL	Ex-Large	Up to 13"				16	



RM Flo-Form Complete Care Boot

Measuring Guide

LCODE: L4396/L4397

Item Code	Size	Foot Length	Infant Shoe	Kids Shoe	Youth Shoe	Men's Shoe	Women's Shoe
20747-S	Small	Up to 8"			1.5		4
20747-M	Medium	Up to 9 1/2"				6.5	8
20747-L	Large	Up to 11"				11	12
20747-XL	Ex-Large	Up to 13"				16	



EZ Stride Anterior Carbon Fiber AFO

Measuring Guide

LCODE: **L1932**

Item Code	Size	Footplate Length	Height Length
CFA-3XSR	XXX-Small Right	7	13"
CFA-3XSL	XXX-Small Left	7	13"
CFA-2XSR	XX-Small Right	8	14"
CFA-2XSL	XX-Small Left	8	14"
CFA-XSR	Ex- Small Right	9"	14
CFA-XSL	Ex- Small Left	9"	14
CFA-SR	Small Right	10"	15
CFA-SL	Small Left	10"	15
CFA-MR	Medium Right	11	16
CFA-ML	Medium Left	11	16
CFA-LR	Large Right	12	17"
CFA-LL	Large Left	12	17"
CFA-Kit	Universal Pad Kit	<i>Suggested code L2820</i>	



* DISTRIBUTED BY RESTORATIVE MEDICAL

EZ Stride Posterior Carbon Fiber AFO

Measuring Guide

LCODE: **L1951**

Item Code	Size	Footplate Length	Height Length
CFP-XSR	Ex- Small Right	9"	14
CFP-XSL	Ex- Small Left	9"	14
CFP-SR	Small Right	10"	15
CFP-SL	Small Left	10"	15
CFP-MR	Medium Right	11	16
CFP-ML	Medium Left	11	16
CFP-LR	Large Right	12	17"
CFP-LL	Large Left	12	17"



* DISTRIBUTED BY RESTORATIVE MEDICAL