


# 2017 PRODUCT GUIDE

HAND		RM Restorative Hand	RM Prosperity Hand	NeuroFlex® Restorative Flex Hand	NeuroFlex® Restorative ThumbEase	RM Cock-Up Wrist Splint	RM Cock-Up Hand Splint with Finger Extension Roll	RM BendEase Hand	RM Pediatric HyperHand
Item #		#20736	#20732	#20721	#20737 / #20740	#20739	#20739cone	#20715	#20717
Hand / Wrist / Forearm	Mild/Moderate Contracture								
	Severe Contracture			Moderate to Closed Palm	Severe Thumb				
Pediatric Sizes Available									
Finger	Severe Flexion in 4th & 5th Digit								
	Severe Thumb ADD								
Moldable By Hand									
Motor or Movement Disorders Related To:									
Cerebral Palsy									
Multiple Sclerosis									
Arthritis									
CVA or Stroke									
Carpel Tunnel Syndrome		Remold Wrist	Remold Wrist	Remold Wrist	Remold Wrist	Molded to 30 degrees	Molded to 30 degrees		

\* This chart is intended to provide caregivers a snapshot of different product features. It may not be comprehensive and should not be used as the sole factor in product selection and use. Each Individual's need are unique. A licensed healthcare practitioner should always verify that the product is appropriate for the individual prior to use.



## RM Restorative Hand Measuring Guide

LCODE: **L3807/3809**

Item Code	Size	Width of MCP	Wrist to Finger
20736XSR	Ex Small -Right	2 3/4" to 3"	Up to 5"
20736XSL	Ex Small -Left	2 3/4" to 3"	Up to 5"
20736SR	Small -Right	Up to 3"	Up to 6 1/2"
20736SL	Small -Left	Up to 3"	Up to 6 1/2"
20736MR	Medium -Right	Up to 3 1/2"	Up to 7"
20736ML	Medium -Left	Up to 3 1/2"	Up to 7"
20736LR	Large -Right	Up to 3 3/4"	Up to 8"
20736LL	Large -Left	Up to 3 3/4"	Up to 8"
20736RZ	Custom - Right		
20736LZ	Custom - Left		



### RM Prosperity Hand

#### Measuring Guide

LCODE: **L3807/3809**

Item Code	Size	Width of MCP	Wrist to Finger
20732SR	Small -Right	Up to 3"	Up to 6 1/2"
20732SL	Small -Left	Up to 3"	Up to 6 1/2"
20732MR	Medium -Right	Up to 3 1/2"	Up to 7"
20732ML	Medium -Left	Up to 3 1/2"	Up to 7"
20732LR	Large -Right	Up to 3 3/4"	Up to 8"
20732LL	Large -Left	Up to 3 3/4"	Up to 8"



### NeuroFlex® Restorative Flex Hand

#### Measuring Guide

LCODE: **L3807/3809**

Item Code	Size	Width of MCP
20720TR	Toddler -Right	2" to 2 1/4"
20720TL	Toddler -Left	2" to 2 1/4"
20720CR	Child-Right	2 1/4" to 2 1/2"
20720CL	Child-Left	2 1/4" to 2 1/2"
20720YR	Youth -Right	2 3/4" to 3"
20720YL	Youth -Left	2 3/4" to 3"



Universal Adult Size LF or RT with cones or cylinders)

20721L	Left w/Regular Cones
20721R	Right w/Regular Cones
20721LRC	Left w/Reverse Cones
20721RRC	Right w/Reverse Cones
20721LCY	Left w/Cylinders
20721RCY	Right w/Cylinders
20721LZ	Custom C-BAR Left
20721RZ	Custom C-BAR Right



### NeuroFlex® Restorative ThumbEase (Pediatric)

#### Measuring Guide

LCODE: L3807/3809

Item Code	Size	Width of MCP	Wrist to Finger	Wrist to Forearm
20740TR	Toddler -Right	2" to 2 1/4"	Up To 3 1/2"	3 1/2"
20740TL	Toddler -Left	2" to 2 1/4"	Up To 3 1/2"	3 1/2"
20740CR	Child-Right	2 1/4" to 2 1/2"	Up To 4"	3 3/4"
20740CL	Child-Left	2 1/4" to 2 1/2"	Up To 4"	3 3/4"
20740YR	Youth -Right	2 3/4" to 3"	Up To 5"	4"
20740YL	Youth -Left	2 3/4" to 3"	Up To 5"	4"



### NeuroFlex® Restorative ThumbEase (Adult)

#### Measuring Guide

LCODE: L3807/3809

Item Code	Size	Width of MCP	Wrist to Finger
20737XSR	Ex Small -Right	2 3/4" to 3"	Up to 5"
20737XSL	Ex Small -Left	2 3/4" to 3"	Up to 5"
20737SR	Small -Right	Up to 3"	Up to 6 1/2"
20737SL	Small -Left	Up to 3"	Up to 6 1/2"
20737MR	Medium -Right	Up to 3 1/2"	Up to 7"
20737ML	Medium -Left	Up to 3 1/2"	Up to 7"
20737LR	Large -Right	Up to 3 3/4"	Up to 8"
20737LL	Large -Left	Up to 3 3/4"	Up to 8"
20737RZ	Custom -Right		
20737LZ	Custom -Left		



### RM Cock-Up Wrist Splint

#### Measuring Guide

Item Code	Size	Width of MCP
20739XSR	Ex Small -Right	2 3/4" to 3"
20739XSL	Ex Small -Left	2 3/4" to 3"
20739SR	Small -Right	Up to 3"
20739SL	Small -Left	Up to 3"
20739MR	Medium -Right	Up to 3 1/2"
20739ML	Medium -Left	Up to 3 1/2"
20739LR	Large -Right	Up to 3 3/4"
20739LL	Large -Left	Up to 3 3/4"



### RM Cock-Up Hand Splint with Finger Extension Roll

#### Measuring Guide

LCODE: L3807/3809

Item Code	Size	Width of MCP
20739rollXSR	Ex Small -Right	2 3/4" to 3"
20739rollXSL	Ex Small -Left	2 3/4" to 3"
20739rollSR	Small -Right	Up to 3"
20739rollSL	Small -Left	Up to 3"
20739rollMR	Medium -Right	Up to 3 1/2"
20739rollML	Medium -Left	Up to 3 1/2"
20739rollLR	Large -Right	Up to 3 3/4"
20739rollLL	Large -Left	Up to 3 3/4"



### RM BendEase Hand

#### Measuring Guide

LCODE: L3807/3809

Item Code	Size	Width of MCP	Wrist to Finger
20715SR	Small -Right	Up to 3"	Up to 6 1/2"
20715SL	Small -Left	Up to 3"	Up to 6 1/2"
20715MR	Medium -Right	Up to 3 1/2"	Up to 7"
20715ML	Medium -Left	Up to 3 1/2"	Up to 7"
20715LR	Large -Right	Up to 3 3/4"	Up to 8"
20715LL	Large -Left	Up to 3 3/4"	Up to 8"



### RM BendEase Hand (Universal, Left or Right)

#### Measuring Guide

LCODE: L3807/3809

Item Code	Size
20715L	Universal Left
20715R	Universal Right



### RM BendEase Hand-Fluidized (Flo-Form™)

#### Measuring Guide

LCODE: L3807/3809

Item Code	Size	Width of MCP	Wrist to Finger
20715SR-F	Small -Right	Up to 3"	Up to 6 1/2"
20715SL-F	Small -Left	Up to 3"	Up to 6 1/2"
20715MR-F	Medium -Right	Up to 3 1/2"	Up to 7"
20715ML-F	Medium -Left	Up to 3 1/2"	Up to 7"
20715LR-F	Large -Right	Up to 3 3/4"	Up to 8"
20715LL-F	Large -Left	Up to 3 3/4"	Up to 8"



### RM Pediatric HyperHand

#### Measuring Guide

LCODE: L3807/3809

Item Code	Size	Width of MCP	Wrist to Finger	Wrist to Forearm
20717IR	Infant -Right	1 3/4" to 2"	Up To 2"	2 1/4"
20717IL	Infant -Left	1 3/4" to 2"	Up To 2"	2 1/4"
20717TJR	Toddler Jr. -Right	1 3/4" to 2"	Up To 3"	2 3/4"
20717TJL	Toddler Jr. -Left	1 3/4" to 2"	Up To 3"	2 3/4"
20717TR	Toddler -Right	2" to 2 1/4"	Up To 3 1/2"	3 1/2"
20717TL	Toddler -Left	2" to 2 1/4"	Up To 3 1/2"	3 1/2"
20717CR	Child-Right	2 1/4" to 2 1/2"	Up To 4"	3 3/4"
20717CL	Child-Left	2 1/4" to 2 1/2"	Up To 4"	3 3/4"
20717YR	Youth -Right	2 3/4" to 3"	Up To 5"	4"
20717YL	Youth -Left	2 3/4" to 3"	Up To 5"	4"

