


EZ Stride Anterior and Posterior Carbon Fiber AFO's

Prevent - Protect - Restore



EZ Stride Anterior Carbon Fiber AFO LCode: 1932

- Lateral strut spirals to the front, providing support for patients who pronate (valgus or eversion)
- One-piece design with “Y” front trim line provides relief to the tibia’s tuberosity area
- Flexible calf strap accommodates a wide range of leg sizes
- Structurally reinforced in high-stress areas
- Optional soft interface kit available (suggested L-code: L2820)

EZ Stride Posterior Carbon Fiber AFO LCode: 1951

- True spiral design
- Lateral strut spirals to medial aspect
- One piece...continuous strand design
- Dynamic footplate
- Super flexible calf band accommodates a wide range of circumferences
- Structurally reinforced in high stress areas

*Distributed by Restorative Medical

EZ Stride Anterior Carbon Fiber AFO	EZ Stride Posterior Carbon Fiber AFO	Size	Footplate Length	Height Length	Men's Shoe	Women's Shoe
CFA-3XSR		XXX-Small Right	7 3/8"	13"		
CFA-3XSL		XXX-Small Left	7 3/8"	13"		
CFA-2XSR		XX-Small Right	8 3/8"	14"		
CFA-2XSL		XX-Small Left	8 3/8"	14"		
CFA-XSR	CFP-XSR	Ex-Small Right	9"	14 1/2"	6.5	8
CFA-XSL	CFP-XSL	Ex- Small Left	9"	14 1/2"	6.5	8
CFA-SR	CFP-SR	Small Right	10"	16"	8	9.5
CFA-SL	CFP-SL	Small Left	10"	16"	8	9.5
CFA-MR	CFP-MR	Medium Right	11 3/8"	16 1/2"	12	13
CFA-ML	CFP-ML	Medium Left	11 3/8"	16 1/2"	12	13
CFA-LR	CFP-LR	Large Right	12 1/4"	17"	15	
CFA-LL	CFP-LL	Large Left	12 1/4"	17"	15	
CFA-Kit		Universal Pad Kit (Suggested code L2820)				



Prevent • Protect • Restore

332 East Broadway • Brandenburg, KY 40108
 (800) 793-5544 • 270-422-5454 • Fax: (270) 422-5453
www.restorativemedical.com

