

# Flo-Form™ Complete Care Boot

## Prevent - Protect - Restore



Boot Size	Foot Size
Small	Up to 8"
Medium	Up to 9.5"
Large	Up to 11"
X-Large	Up to 13"

- Safe and comfortable solution for contractures of the foot and ankle, foot drop, plantar flexion and conditions resulting from various types of neuropathy, inflammation or nerve injury.
- Patented fluidized cushioning creates the perfect blend of support, comfort and contour to help achieve the best possible therapeutic positioning.
- Helps to maintain neutral alignment of the foot in a non-ambulatory patient.
- Fluidized sections to help redistribute force include lower calf and Achilles area.
- Suspends the heel to prevent excessive pressure and skin breakdown.
- Anti-skid pad for safe brief standing and pivot transfers (not designed for ambulation).
- High temperature plastic is molded to help float the heel, and is easily remolded if needed to accommodate moderate to severe loss of motion.
- Soft, durable fabric next to patient's skin provides protection and comfort.
- Anti-rotation bar helps to maintain alignment of lower leg and correct hip rotation.
- Optional adjustable straps maximize dorsiflexion to help prevent foot drop.