



## Measuring Guide

NeuroFlex® Restorative Knee Item Code	NeuroFlex® Restorative HyperExtension Knee Item Code	Size	TTK	KTC	Thigh Cir.	Calf Cir.
20725CS	20725CSH	Child/Slim	5"	6"	7 ½-9 ½"	6 ½-8 ½"
20725C	20725CH	Child	5"	6"	8-12"	7-10"
20725CF	20725CFH	Child/Full	5"	6"	8-13"	7 ½-11"
20725XSS	20725XSSH	Ex-Small/Slim	6"	6"	8-12"	7-10"
20725XS	20725XSH	Ex-Small	6"	6"	8-13"	7 ½-11"
20725XSF	20725XSFH	Ex-Small/Full	6"	6"	10-15"	8-12"
20725SS	20725SSH	Small/Slim	6 ½"	6 ½"	8-13"	7 ½-11"
20725S	20725SH	Small	6 ½"	6 ½"	10-15"	8-12"
20725SF	20725SFH	Small/Full	6 ½"	6 ½"	13-17"	10-15"
20725MS	20725MSH	Medium/Slim	7 ½"	7 ½"	10-15"	8-12"
20725M	20725MH	Medium	7 ½"	7 ½"	13-17"	10-15"
20725MF	20725MFH	Medium/Full	7 ½"	7 ½"	16-22"	14-18"
20725LS	20725LSH	Large/Slim	8 ½"	8 ½"	13-17"	10-15"
20725L	20725LH	Large	8 ½"	8 ½"	16-22"	14-18"
20725Z	20725ZH	Custom				