

RESTORATIVE MEDICAL PRODUCT SUMMARY

All Restorative Medical splints can be adjusted as the patient improves by bending, heat molding, or adjusting hinges.

HAND



CE **RM Restorative Hand**
Product # 20736
L-Code L3807/L3809
Comfortably addresses wrist/hand/finger lost range and tension.



CE **RM BendEase Hand - Universal**
Product # 20715
L-Code L3807/L3809
Moldable aluminum base, used for moderate lost range (universal size fits most adults).



BendEase Hand with Flo-Form™ - Universal
Product # 20715UNFLO
L-Code L3807/L3809
Flo-Form™ fluidized medium provides even pressure distribution at point of orthotic contact.



CE **RM Cock-Up Wrist Splint**
Product # 20739
30 degree wrist cock-up Provides moderate, comfortable support of hand and wrist, post injury or surgery.



CE **NeuroFlex® Restorative ThumbEase**
Product # 20737
L-Code L3807/L3809
For severely adducted thumbs and thenar eminence collapse.



CE **RM BendEase Hand - Sized**
Product # 20715
L-Code L3807/L3809
Bend-by-hand moldable aluminum base, used for moderate lost range.



BendEase Hand with Flo-Form™ - Sized
Product # 20715
L-Code L3807/L3809
Flo-Form™ fluidized medium provides even pressure distribution at point of orthotic contact.



CE **RM Cock-Up Hand Splint with Finger Extension Roll**
Product # 20739roll
L-Code L3807/L3809
30 degree wrist cock-up with removable large roll when finger extension is desired.

FOOT/ANKLE




CE **Flo-Form™ Ankle Foot Orthosis**
Product # SUN TAFO
L-Code L4396/4397
Helps to maintain neutral alignment of the foot and to relieve pressure on the heel using fluidized technology.



CE **NeuroFlex® Restorative FootEase**
Product # 20742
L-Code L4396/4397
Addresses mild to moderate plantar flexion contractures, extensor tone and plantar fasciitis.



CE **NeuroFlex® Restorative Flex Hand**
Product # 20721
L-Code L3807/L3809
Addresses tightest closed fist to moderate finger extension with graduated size cones, reverse cones or cylinders.



CE **RM Prosperity Hand**
Product # 20732
L-Code L3807/L3809
Slit between the 3rd & 4th digits allows for varying degrees of range of motion of fingers.



Flo-Form™ Complete Care Boot
Product # 20747
L-Code L4396/4397
Helps to maintain neutral alignment of the foot and to relieve pressure on the heel using fluidized technology.



CE **RM Safe Boot I**
Product # 207441
L-Code L4398
Addresses flaccid foot drop and suspends the heel.



EZ Stride Posterior Carbon Fiber AFO
Product # CFP
L-Code L1951
Enhance normal gait for persons with flaccid foot drop. The lateral strut design spirals to medial aspect. Best for patients who cannot tolerate tibial pressure.
*Distributed by Restorative Medical



EZ Stride Anterior Carbon Fiber AFO
Product # CFA
L-Code L1932
Enhance normal gait for persons with flaccid foot drop. The lateral strut design spirals to anterior aspect providing support for patients who pronate (valgus or eversion). Tibial pad available.
*Distributed by Restorative Medical



CE **NeuroFlex® Restorative Flex Boot**
Product # 20734
L-Code L4396/4397
Addresses mild to severe lost range and abnormal muscle tone.



CE **RM Safe Boot II**
Product # 20744II
L-Code L4396/L4397
Addresses plantar flexion, flaccid foot drop, plantar fasciitis and mild abnormal muscle tone.



RESTORATIVE MEDICAL PRODUCT SUMMARY

The most comprehensive selection of static splints with dynamic properties and technologies that provide therapeutic positioning and aid in the management and rehabilitation of contractures in sizes from infants to adults.



Prevent • Protect • Restore

restorativemedical.com

BACK/NECK



CE **RM VertebrEase TLSO**
Product # 20730
L-Code L0456/L0457

Helps to improve posture including kyphosis, lordosis, scoliosis and other spinal deformities. Padded shoulder straps encourage upright posture.



CE **NeuroFlex® Restorative Safe Spine**
Product # 20728
L-Code L0454/L0455

Comfortable spine support for conditions like kyphosis, scoliosis, lordosis. Padded shoulder straps encourage upright posture.



CE **NeuroFlex® Restorative Cervical HyperExtension Stop**
Product # 20728H
L-Code L0999

Fits into our TLSO and addresses neck hyper-extension, head leaning.



CE **RM VertebrEase LSO**
Product # 20729
L-Code L0631/L0648

Helps to immobilize and support the lumbar spine for soft tissue injuries and relief of low back pain. Shown with optional shoulder straps. Available with pendulous abdomen front.



CE **NeuroFlex® Restorative Cervical Extension, Forward**
Product # 20728C
L-Code L0999

Fits into our TLSO and addresses forward head upper body posture.



CE **NeuroFlex® Restorative Safe Cerv, Left or Right**
Product # 20728SCL/R
L-Code L0999

Fits into our TLSO and addresses left or right head leaning with or without forward leaning.

KNEE/HIP



CE **NeuroFlex® Restorative Knee**
Product # 20725
L-Code L1831

Addresses mild to severe flexion contractures, including those with below knee amputations (BKA).



Knee Separator (Accessory to NeuroFlex® Restorative Knee)
Product # 20725SP

Helps to gradually separate the knee joint before transitioning to the NeuroFlex® Restorative Knee.



CE **NeuroFlex® Restorative Kentucky Kollar**
Product # 20735
L-Code L0113

Addresses dystonia, torticollis and wry neck conditions.



CE **RM RestAir Knee**
Product # 20755air
L-Code L1831/L1847

Addresses knee flexion contractures with air technology.



CE **NeuroFlex® Restorative Hyperextension Knee**
Product # 20725H
Suggested L-Code L1831

Helps to improve and/or restore extension lost range of the knee with effective prolonged low load passive stretch.



CE **RM RestAir Hip**
Product # 20760air
L-Code L1652

Comfortably and safely helps to address hip adduction & scissoring.

ELBOW



CE **Neuroflex® Restorative Elbow**
Product # 20724
L-Code L3760

Addresses mild to severe flexion contractures, including those with below elbow amputations (BEA).



Elbow Separator (Accessory for NeuroFlex® Restorative Elbow)
Product # 20724SP

Helps to gradually separate the elbow joint before transitioning to the NeuroFlex® Restorative Elbow.



CE **RM RestAir Elbow**
Product # 20724air
Suggested L-Code L3760

Addresses elbow flexion contractures with air technology.



CE **NeuroFlex® Restorative Hyper-Extension Elbow**
Product # 20724H
L-Code L3760

Helps to improve and/or restore extension lost range of the elbow