


HAND



Measuring Guide

RM BendEase Hand - Sized Item Code	RM BendEase Hand - Universal Item Code	Size	Width of MCP	Wrist to Finger
20715SR		Small-Right	Up to 3"	Up to 6 ½"
20715SL		Small-Left	Up to 3"	Up to 6 ½"
20715MR		Medium-Right	Up to 3 ½"	Up to 7"
20715ML		Medium-Left	Up to 3 ½"	Up to 7"
20715LR		Large-Right	Up to 3 ¾"	Up to 8"
20715LL		Large-Left	Up to 3 ¾"	Up to 8"
	20715L	Universal-Left		
	20715R	Universal-Right		