



## Measuring Guide

Restorative Hand Item Code	NeuroFlex® Restorative ThumbEase Item Code	Size	Width of MCP	Wrist to Finger
20736XSR	20737XSR	Ex Small-Right	2 ¾" to 3"	Up to 5"
20736XSL	20737XSL	Ex Small-Left	2 ¾" to 3"	Up to 5"
20736SR	20737SR	Small-Right	Up to 3"	Up to 6 ½"
20736SL	20737SL	Small-Left	Up to 3"	Up to 6 ½"
20736MR	20737MR	Medium-Right	Up to 3 ½"	Up to 7"
20736ML	20737ML	Medium-Left	Up to 3 ½"	Up to 7"
20736LR	20737LR	Large-Right	Up to 3 ¾"	Up to 8"
20736LL	20737LL	Large-Left	Up to 3 ¾"	Up to 8"
20736XLR	20737XLR	Ex Large-Right	Up to 4"	Up to 9"
20736XLL	20737XLL	Ex Large-Left	Up to 4"	Up to 9"
20736RZ	20737RZ	Custom-Right		
20736LZ	20737LZ	Custom-Left		